




Impression taking

	Article	Art. No.	Ø
 <p>3 mm 10 mm</p>	<p>CONELOG® Impression post, open tray incl. fixing screw (the fixing screw can be shortened extra-oral by 3 mm with a screwdriver, hex)</p> <p>Material Titanium alloy</p>	C2121.3300	3.3 mm
		C2121.3800	3.8 mm
		C2121.4300	4.3 mm
		C2121.5000	5.0 mm
 <p>10.7 mm</p>	<p>CONELOG® Impression post, closed tray incl. impression cap, bite registration cap and fixing screw</p> <p>Material Titanium alloy / POM</p>	C2110.3300	3.3 mm
		C2110.3800	3.8 mm
		C2110.4300	4.3 mm
		C2110.5000	5.0 mm
	<p>Impression cap for impression post, closed tray (5 units)</p> <p>Material POM</p>	J2111.3300	3.3 mm
		J2111.3800	3.8 mm
		J2111.4300	4.3 mm
		J2111.5000	5.0 mm