


SCREW-LINE


Implants with snap-in insertion posts

	Article	Art. No.	Ø	L	A Ø
	CONELOG® SCREW-LINE Implant, Promote® plus incl. snap-in insertion post and cover screw, sterile Material Titanium Grade 4	C1066.3309	3.3 mm	9 mm	2.7 mm
		C1066.3311		11 mm	
		C1066.3313		13 mm	
		C1066.3316		16 mm	
		C1066.3807	3.8 mm	7 mm	3.5 mm
		C1066.3809		9 mm	
		C1066.3811		11 mm	
		C1066.3813		13 mm	
		C1066.3816	16 mm	3.9 mm	
		C1066.4307	7 mm		
		C1066.4309	9 mm		
		C1066.4311	11 mm		
		C1066.4313	13 mm	4.6 mm	
		C1066.4316	16 mm		
		C1066.5007	7 mm		
		C1066.5009	9 mm		
		C1066.5011	11 mm	5.0 mm	
		C1066.5013	13 mm		
C1066.5016	16 mm				

Note

CONELOG® SCREW-LINE implants, Promote® plus, with Art. No. C1066.xxxx and C1065-xxxx can be used exclusively with the insertion tools (Art. No. J5300.0031, J5300.0032, J5300.0033, J5300.0034 or J5300.0035, J5300.0036 and J5300.0037).

Implants with screw-mounted insertion posts

	Article	Art. No.	Ø	L	A Ø
	CONELOG® SCREW-LINE Implant, Promote® plus incl. screw-mounted insertion post and cover screw, sterile Material Titanium Grade 4	C1065.3309	3.3 mm	9 mm	2.7 mm
		C1065.3311		11 mm	
		C1065.3313		13 mm	
		C1065.3316		16 mm	
		C1065.3807	3.8 mm	7 mm	3.5 mm
		C1065.3809		9 mm	
		C1065.3811		11 mm	
		C1065.3813		13 mm	
		C1065.3816	16 mm	3.9 mm	
		C1065.4307	7 mm		
		C1065.4309	9 mm		
		C1065.4311	11 mm		
		C1065.4313	13 mm	4.6 mm	
		C1065.4316	16 mm		
		C1065.5007	7 mm		
		C1065.5009	9 mm		
		C1065.5011	5.0 mm	11 mm	4.6 mm
		C1065.5013		13 mm	
		C1065.5016		16 mm	

Note

Implants with the screw-mounted insertion post (Art. No. C10655.xxxx) are to be used for template-guided implant insertion with the SCREW-LINE Guide System.

Note

CONELOG® SCREW-LINE implants, Promote® plus, with Art. No. C1066.xxxx and C1065-xxxx can be used exclusively with the insertion tools (Art. No. J5300.0031, J5300.0032, J5300.0033, J5300.0034 or J5300.0035, J5300.0036 and J5300.0037).