



Scanbodies

	Article	Art. No.	Ø
	<p>CONELOG® Scanbody* for optical, 3-dimensional localization of CONELOG® Implants in the mouth or CONELOG® Lab analogs in the working model, incl. abutment screw, sterile</p> <p>Not compatible with the CEREC and inLab systems from Dentsply Sirona</p> <p>Material PEEK</p>	C2600.3310	3.3 mm
		C2600.4310	3.8 mm
			4.3 mm
		C2600.5010	5.0 mm
	<p>CONELOG® ScanPost for Sirona® Scanbody for digital recording of the CONELOG® Implant or lab analog position and for further processing in the CEREC and inLab systems from Dentsply Sirona, incl. abutment screw</p> <p>Material Titanium alloy</p>	C2620.3306	3.3 mm
		C2620.3806	3.8 mm
		C2620.4306	4.3 mm
		C2620.5006	5.0 mm

* Please check whether the CONELOG® Scanbody is available in the CAD software used.
CAD libraries for selected CONELOG® Prosthetic components are available for free download here:
www.camlog.com/en/media-center/cad-libraries

Matching Sirona® Scanbodies size S for CONELOG® ScanPost and CONELOG® Titanium base CAD/CAM, crown, with Ø 3.3/3.8/4.3 mm:

For Omnicam®: Article number 6431311
For Bluecam®: Article number 6431295

Matching Sirona® Scanbodies size L for CONELOG® ScanPost and CONELOG® Titanium base CAD/CAM, crown, with Ø 5.0 mm:

For Omnicam®: Article number 6431329
For Bluecam®: Article number 6431303

Sirona® Scanbodies are available from Dentsply Sirona.